



QUALITY ASSURANCE ISO 9001:2015



ASV Mphi VET

Portable. Ergonomic. Easy to use.





Mphi **VET**

LASER THERAPY



Device Specifications:

Pre-settled Programs	18 programs for canine and feline 21 Programs for equine
MLS Modulation Modes	Continuous Pulsed Waves (CP5) Frequenced Pulsed Waves (FPS)
Modulation Frequency	1 to 2000 Hz with 1 Hz steps
Power Levels	25%, 50%, 75% and 100%
Fixed Duty Cycle	50% or Variable
Treatment Duration	1" to 99'59" with 1" Steps
Handpiece	With optictal group composed of 1 MLS Source Power up to 1.1 W (Peak Power: 25 W)
Target Area	Ø 2 cm Highlighted by high-efficiency red light LEDs
Display	High resolution LCD color touchscreen user interface
Other Features	Simultaneous management of 2 independent channels Possibility to customize every treatment parameter Personal area to save custom protocols Automatic calculation of emitted energy according to the settled parameters